INSTRUCTIONS FOR USE

ASAHI REXEED™-A SERIES DIALYZER

1. INDICATION FOR USE

Asahi REXEED™-A series dialyzer (hereafter called REXEED-A) is intended for use in hemodialysis (HD), hemodiafiltration (HDF) and hemofiltration (HF) for the treatment of patients suffering from acute or chronic renal failure.

2. CONTRAINDICATIONS

None known. In general, contraindications for HD, HDF and HF are applicable.

3. ADVERSE REACTIONS

- Consult a responsible physician immediately whenever an adverse reaction occurs.
- Some patients may experience hypersensitivity reactions during treatment. Symptoms and signs include acute shortness of breath with wheezing; respiratory arrest; itching; flushing; mood discomfort; abnormal sweating; hives; generalized or localized skin redness; edema of the face, hands, or feet; hypertension; hypotension; elevated pulse rates; arrhythmia; ocular hyperemia; hypoesthesia; fever; leukopenia; and thrombocytopenia.
- Side effects such as hypertension, hypotension, headache, chest pain, stomachache, fatigue and nausea which may be associated with hypovolemia or hypervolemia can usually be avoided by careful management of the patient's fluid volume and electrolyte balance as well as the blood flow rate and transmembrane pressure (TMP).

4. WARNINGS

Follow these instructions and those of the HD machine supplier. Improper use could cause injury or death. The manufacturer will not be responsible for patient safety if these instructions and the HD machine supplier instructions are not followed. Make sure these instructions are available to the operator at the treatment site.

- REXEED-A is intended for single use only. Reprocessing of REXEED-A may lead to adverse reactions for the patient and/or device failure.
- The user must use REXEED-A before the expiration date. Refer to the REXEED-A product label for the expiration date.
- REXEED-A must be kept in a clean, dry environment away from direct sunlight, heat, and humidity and at temperatures of 0-30°C (32-86°F).
- REXEED-A, when shipped by the manufacturer, is filled with fluid. This fluid must not be allowed to freeze. REXEED-A with frozen fluid must never be used, even after the thawing of the fluid.
- REXEED-A is a plastic product, and must not be exposed to vibration or impact during transportation, storage or handling.
- REXEED-A must not be exposed to organic solvents such as alcohol to prevent cracking or deformation.
- Inspect each REXEED-A before use, and replace it if it shows signs of damage, damage or improper sealing of the sterilization bag, leakage of fluid from REXEED-A into the sterilization bag, or misaligned or separated stoppers.
- Ensure that the dialysate and blood compartments of REXEED-A are rinsed before use according to the priming procedure described below.
- Use aseptic technique throughout priming and treatment.
- Do not remove the fluid before rinsing and filling the blood compartment with saline to prevent any air from entering REXEED-A.
- The saline must be a physiological solution that is sterile and non-pyrogenic.
- Due to its high ultrafiltration coefficient, REXEED-A must only be used with HD machines equipped with a UF controller (Refer to the instructions of the HD machine).
- Dialysate must meet the regional or ISO water quality standards for dialysis to minimize the risk of adverse reactions.
- The use of blood and dialysate port connectors (female) attached to blood and dialysate lines that are designed in compliance with ISO 8637:2004 and ISO 8638:2004 is recommended.
- Ensure that the sensors, detectors and alarm system of the HD machine are in proper operating condition including the venous air sensor and the blood leak detector in accordance with the instructions of the HD machine before starting treatment.
- Ensure that air is completely purged from REXEED-A and blood lines during priming; failure in purging air may result in blood clotting during treatment. Further, ensure that no air enters the patient's blood vessel during blood return.
- During priming and treatment, verify that there is no leakage from the dialyzer, blood lines, blood ports, and their connecting parts. If leakage is observed, replace the dialyzer or the blood lines.
- If dialysate is to be used instead of saline or substitution fluid, care must be taken to ensure that the purity of the dialysate is in accordance with the applicable water quality standards (for example, dialysate intended for infusion defined in ISO 11663:2009).
- In the event of a problem during treatment, such as blood leakage or coagulation, immediately discontinue the treatment under the direction of a physician and replace REXEED-A with a new primed dialyzer.
- REXEED-A has multiple connection points. Failure to make the

proper connections can result in compromised treatment, injury, or death due to blood or therapy fluid leak.

- REXEED-A may remove medications given in the arterial blood line. Give medications post dialyzer unless instructed otherwise.
- Do not permit the TMP to exceed 500 mmHg (66 kPa) with TMP being defined as per ISO 8637:2004.
- The recommended ranges of the treatment parameters specified in the APPENDIX should be considered when planning treatment for patients.
- Dispose the used REXEED-A and the blood lines in strict accordance with official and institutional standards for medical waste disposal.

5. SPECIFICATIONS & PERFORMANCE

Membrane, REXBRANE™ (Polysulfone / PVP); Housing and headers, Styrene-butadiene block copolymer; Stoppers, Polyethylene and hydrogenated styrene-butadiene block copolymer; Potting material, Polyurethane

Sterilized by gamma rays

Blood pathway is non-pyrogenic.

Refer to the APPENDIX for design specifications and *in vitro* performance information regarding individual REXEED-A models. In a clinical setting, the performance of REXEED-A may be different from that *in vitro*; this performance may vary with the conditions and duration of treatment.

6. ANTICOAGULATION

Anticoagulant requirements will vary with the patient's condition; the type, quantity, and method of anticoagulant administration should be in accordance with the directions by a physician who is familiar with the patient's condition.

7. OPERATION**1. PRIMING**

The order of procedures 2) and 3) can be reversed.

- 1) Place REXEED-A vertically in the holder, with the red header (arterial end) down.
- 2) Connect the arterial blood line to the saline bag for infusion and fill the arterial blood line completely with saline. Connect the arterial and venous blood lines to the blood inlet and outlet ports. Be careful not to introduce air into the hollow fibers and the blood lines. Rinse the hollow fibers and the blood lines with at least 1000 mL saline at approx. 100 mL/min. If any air is trapped in the dialyzer or blood lines, introduce additional saline for rinsing while gently tapping the dialyzer with one's hands or pinching and releasing the blood lines until the air is completely removed. Do not attempt to remove the air by tapping the dialyzer with hard instruments, such as forceps, because this may damage the dialyzer.
- 3) Connect the dialysate lines to the dialysate ports. Allow the dialysate to flow at approx. 500 mL/min for 2 min or more. For HDF, connect the substitution line to the substitution port and prime the line.
 - For HF, different priming procedure can be obtained on request.
 - Instead of saline, dialysate intended for infusion may be used for priming and ending the treatment.

2. STARTING TREATMENT

- 1) Connect the arterial blood line to the patient's arterial cannula. Start blood pump (approx. 100 mL/min) to displace saline.
- 2) Connect the venous blood line to the venous cannula.
 - For HDF, the ultrafiltration flow rate should be monitored carefully under the direction of a physician who is familiar with the patient's condition. In some patients, for example, those with a high hematocrit value, a high ultrafiltration flow rate may cause elevation of the blood inlet pressure.

3. ENDING TREATMENT

- 1) Adjust the UF flow rate and stop the blood pump.
- 2) Remove the arterial cannula from the patient and connect the arterial blood line to the saline bag. Feed 100-200 mL saline into the hollow fibers and the blood lines by starting the blood pump slowly to rinse back patient's blood.
- 3) Clamp the venous blood line and remove the venous cannula from the patient.

8. WARRANTY AND LIMITED LIABILITY

- The manufacturer warrants that REXEED-A is manufactured in accordance with its specifications and in compliance with ISO 9001:2008, ISO 13485:2003, applicable industry standards, and regulatory requirements.
- If any defects caused during manufacturing or transportation are brought to the attention of the manufacturer's representative, the manufacturer will replace REXEED-A with a new one through the representative.
- The manufacturer will not be liable for any damage caused by the reuse of REXEED-A, misuse, improper handling, non-compliance with warnings and instructions, damage caused by events occurring after the product is released, failure to ensure the product is in proper condition before use, or any warranty given by independent distributors.

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ASAHI KASEI MEDICAL CO., LTD.
ASAHI KASEI MEDICAL EUROPE GmbH