


Technical Data; BK-F series Filtryzer


Type	BK-1.3F	BK-1.6F	BK-2.1F
Housing	Material: Polystyrene		
Fibers	Material: Polymethylmethacrylate		
	Inner diameter (μm): 200		
	Membrane thickness (μm): 30		
	Effective surface area (m ²): 1.3	1.6	2.1
Potting Material	Polyurethane		
Sterilization	Gamma-ray Irradiation		
Blood Volume (mL)	76	94	126
Clearance <i>in vitro</i> (mL/min)*			
Urea	designed 183	190	195
Creatinine	designed 160	172	181
Phosphate	designed 144	157	171
Vitamin B ₁₂	designed 99	111	128
Inulin	designed 53	61	73
UFR <i>in vitro</i> (mL/hr, at 13.3kPa (100mmHg))**	1,600	2,000	2,600
Max. TMP (kPa (mmHg))	66 (500)		

* Clearances are data with aqueous solution
 Q_s: 200 ±4mL/min, Q_o: 500 ±10mL/min, Q_r: 10 ±2 mL/min, Temp.: 37 ±1°C
 ** UFRs are typical data with bovine blood. (Ht 30 ±3%, TP 6 ±0.5g/dL)
 Q_s: 200 ±4mL/min, TMP: 13.3 ±1.3kPa (100 ±10mmHg), Temp.: 37 ±1°C

Clearance



UFR



TORAY

CE 0123

EC REP

Toray International Italy S.r.l.
 Via Mecenate 86, 20138 Milan, ITALY
 TEL: 39-02-580-39133 / FAX: 39-02-580-16317

Exporter:
Toray Medical Co., Ltd.
 Dialysis Products Business Division
 8-1, Mihama 1-chome, Urayasu, Chiba 279-8555, JAPAN
 TEL: 81-47-700-7537 / FAX: 81-47-700-7558 / E-MAIL: TMC_INTL_FL@tmc.toray.co.jp

Manufacturer:
Toray Industries, Inc.
 1-1, Nihonbashi-Muromachi 2-chome, Chuo-ku, Tokyo 103-8666, JAPAN
 TEL: 81-3-3245-5144 / FAX: 81-3-3245-5609

TORAY

FILTRYZER® BK-F SERIES

Hollow Fiber Dialyzer



PMMA for better quality of life


The PMMA membrane offers excellent clinical benefits to renal failure patients.




PMMA Membrane; More Efficient for the Elimination of Large Molecular Weight Substances.

It is suggested that pathogenic substances, such as erythropoiesis inhibitory factors^{1-4,7)}, substances related to bone disease⁸⁾, immunosuppressive factors⁹⁾ etc, exist in the large molecular weight range from 50,000 to 1,000,000 daltons. The large molecular weight substances, however, are not removed in conventional HD treatment so as not to lose excessive albumin.

With that in mind, Toray developed the BK-F series to have a larger pore size than the conventional BK-U and -P series and it is expected to remove larger molecular weight substances than albumin.



BK-F Removes Large Molecular Weight Substances⁷⁾




Concentrated dialysate of chronic renal failure (CRF) patients who were dialyzed with the new BK-F membrane was fractionated with a Sephacryl gel column.


Substances in the concentrated filtrate were separated by molecular weight into three fractions called KR4-0, KR4-1, and KR4-2.

Fractions of KR4-1 and 4-2 primarily consist of IgG and albumin, respectively. KR4-0 fraction whose estimated molecular weight was between 500,000 and 1,000,000 daltons was found mainly in the dialysate when BK-F membrane was employed.

An Inhibitor of Erythropoiesis Exists in Large Molecular Weight Substances?



Hematocrit Change Using BK-F Membrane⁷⁾



A clinical study using BK-F hemodialyzers was performed on 8 CRF patients without administering erythropoietin (EPO). As illustrated above, 6 out of 8 CRF patients maintained hematocrit (Hct) at an acceptable level (25-30%) for 11 months, and 2 out of 8 patients, who had extremely low Hct before the study, showed improvement in Hct from 21-22% to 25% after 11 months.

Minimized Loss of Total Protein⁷⁾

	Control	Duration of Dialysis using BK-F	
		3-5 Months	8-10 Months
Albumin	3.8 ± 0.2 (16)	3.7 ± 0.3 (24)	3.7 ± 0.3 (24)
IgG	1340 ± 347 (8)	1393 ± 385 (8)	1325 ± 360 (8)
IgA	114.5 ± 84.2 (8)	135.9 ± 96.7 (8)	120.9 ± 84.3 (8)
IgM	102.3 ± 50.7 (8)	111.0 ± 57.0 (8)	100.5 ± 47.7 (8)
Total Cholesterol	159.2 ± 25.5 (8)	190.8 ± 31.5 (24)	181.2 ± 31.5 (16)

NS: Not Significant, p<0.05: Significant difference.

Even though large molecular weight substances can be removed, the sieving coefficient for albumin is only about 0.03.


As illustrated above (mean ±SD (n)), the elution of albumin, IgG, IgA, and IgM were studied. No significant changes were obtained during 10 months using BK-F in the HD mode.

Long-term Multicenter Clinical Evaluation

BK-F dialyzers remove β_2 -microglobulin, which is considered to be one of the pathogenic substances of amyloidosis, as our conventional BK-U and -P series can.

It has been reported in the long-term multicenter studies both at home and abroad with the BK series that it is effective in preventing carpal tunnel syndrome (CTS) in long-term dialysis patients and in improving itching¹¹⁾, joint pain, bone cysts¹⁰⁾, etc.

Time Course of Changes of Pain Index during BK Membrane Hemodialysis¹¹⁾



REFERENCES:

- Dapner TA, et al. Effective of Low Dose Erythropoietin: A Possible Advantage of High Flux Hemodialysis. *Trans Am Soc Artif Intern Organs* 1990; 36: M223.
- Ueno K, Watanabe E, Aoyagi H, et al. Kidney and Dialysis 1990; 28 (Suppl) 90-93 (in Japanese).
- Ohmura K, Sawaoka K, Chiba E, et al. Kidney and Dialysis 1991; 30 (Suppl) 128-132 (in Japanese).
- Sakashita K, Tsutsui T, Yamamoto N, et al. Influence of Hemodialysis Membranes on Clinical Effect of Recombinant Human Erythropoietin. *Jpn.J.Artif.Organs* 1992; 21: 855-860.
- Shimizu M, Kumegawa M. Japanese Patent Laid Open 1989-287099.
- Hörl WH, et al. Physicochemical characterization of a polypeptide present in uremic serum that inhibits the biological activity of polymorphonuclear cells; *Proc. Natl. Acad. Sci* 1990; 87: 6353-6357.
- Kobayashi H, Ono T, Yamamoto M, et al. Removal and characterization of hemopoiesis inhibitors using a large-pore membrane. *Kidney and Dialysis* 1993; 34 (Suppl) 154-157 (in Japanese).
- Hakim RM, et al. *Kidney International* 1984; 26: 194-200.
- Tetta C, Camussi G, Turello E, et al. Production of

INSTRUCTIONS:

Filtrizer BK-F series is a medical device intended for hemodialysis (HD), but must not be used for HDF (hemodiafiltration) or HF (hemofiltration) due to the higher permeability of larger molecular weight proteins such as albumin.

This device must be used by or at the direction of a physician.

Patients with bleeding tendencies or coagulation disorders must be carefully evaluated by the physician. When adverse reactions are observed, the patients must

be promptly treated under the direction of the physician. For some reactions, manipulation of blood flow rate, ultrafiltration rate, and electrolytic balance can be applied.

The "Instructions for Use" should be read thoroughly prior to the use of this medical device.

Filtrizer is manufactured in accordance with "Approval Standard of Artificial Kidney" by the Ministry of Health, Labour and Welfare of Japanese Government.

Each unit is carefully tested, sterilized and packaged prior to shipment. Toray cannot assume any responsibility


Biocompatibility

Polymethylmethacrylate (PMMA) membranes are less complement activating, have less reduction in neutrophil counts, and do little damage to platelets¹²⁾ during dialysis, because of its superior biocompatibility.

In clinical use of PMMA membrane, the plasma levels of cytokines, such as TNF- α and IL-1 β were consistently reduced in respect to the predialytic values relative to other membranes⁹⁾.

Thanks to gamma-ray sterilization, there is no concern about side effects due to residual sterilants.

Change of C3a During Dialysis⁹⁾



Cytokines in Hemodialysis. *Blood Purif.* 1990; 8: 337-346.

10) Nakamoto M, Goya T, Takahashi H, et al. The 37th Congress of Japanese Society for Dialysis Therapy 1992: 288.

11) Arakawa M, et al. Long-Term Multicenter Study on β_2 -Microglobulin Removal by PMMA BK Membrane. *Nephrol Dial Transplant* 1991; 6 (Suppl 2) 69-74.

12) Akizawa T, Nishiyama H, Koshikawa S. Plasma β_2 -microglobulin levels in chronic renal failure patients. *Int. Soc. Art. Organs* 1981; 5: 54-58.

for damage that may occur during transport or due to mishandling.

Filtrizer is filled with sterile water. Before starting dialysis, rinse it out with one liter or more of physiological saline solution.

Filtrizer is designed for single use only.

Since Filtrizer BK-F series has high ultrafiltration rates, it is necessary to use a dialysis machine equipped with a volumetric ultrafiltration rate controller.