

## ComfortLine Canada: ...special and functional

The ComfortLine Canada therapy chair is especially designed to offer higher comfort for your patient. The unique headrest with tiltable sides and neck pillow allow a comfortable and safe head position. The specially shaped back and sitting segments guarantee a convenient and long-lasting seating position. Besides the ComfortLine Canada offers all advantages of the well known standard Comfort-Line, for example the 5 easy to access motors, the removable armrest and the activation of the shock position by foot.

Different colours make the chair fit perfectly into the desired environment.



Standard colour:

Chrom grey 

Special colours for  
bi-colour design:

Feuer 

Atoll 

Lima 

Amazonas 

## Technical Data

### Dimensions

Seat section width	60 cm
Chair width incl. arm rests	90 cm
Total length	215 cm
Total weight	92 kg

### Adjustments

Back section	75° - 0°
Seat section	30° - 5°
Leg section	5° - 30°
Foot rest electrical	25 cm
Foot rest manual	18 cm
Height adjustment	63 - 83 cm

### Safety

Safe maximum load	180 kg
Angle of shock position	- 12°
Time for shock position	10 sec

### Electrical specifications

Power supply / voltage	120 V 50/60 Hz
Voltage	173 VA
Protection class	2
Protection category	IPX 4
In compliance with	EN 60601
Motors	5



Dialysis  
Transfusion  
Oncology



Right for technical modifications reserved  
B-MOBCLCANIE01



Bionic Medizintechnik GmbH ... a JMS-corporation • Max-Planck-Straße 21  
D-61381 Friedrichsdorf/Germany • Phone +49 (0) 6172-7576-0  
Fax +49 (0) 6172-7576-10 • info@bionic-jms.com • www.bionic-jms.com  
Chief Medical Supplies Ltd. - 411-19th Street S.E, Calgary, Alberta T2E 6J7  
Phone: (866) 620-6034 • Fax: (403) 207-6074 • contact@chiefmedical.com  
www.chiefmedical.com





The ergonomically shaped upholstery adapts in an optimal way to the patient's body and the water resistant surfaces allow easy cleaning.

Five motors make the chair adjustable into all positions desired by the patient and by the nursing staff. Removable armrests and an adjustable foot and head section make the long term dialysis treatment more comfortable for the patient.

The ComfortLine provides very high stability whilst in horizontal position. The powerful motor system allows the shifting to shock position at high speed. The activation of the shock position by foot pedal ensures that immediate measures can be initiated without pressing the emergency button by hand first.

A height adjustable pillow provides optimal head support. The telescopic footrest can be adjusted electrically to give very small up to very tall persons a secure feeling during treatment.

ComfortLine is very well known in Europe since many years. For Canada we have developed a unique shaped upholstery. All other known features of the ComfortLine remain identical. For further information please visit our website:

[www.therapychair.com](http://www.therapychair.com)