



## Chief Medical Supplies Ltd.

411 - 19 Street S. E.  
Calgary, Alberta, Canada  
T2E 6J7

Telephone: +1 403.207.6034 Facsimile: +1 403.207.6074

### CITRASATE FORMULATION LIST

*(in mmole/ Litre after standard 1:45 dilution)*

Formulation	Sodium	Potassium	Calcium	Magnesium	Chlorides	Dextrose	Citric Acid	Sodium Acetate	DIN
C0125	100.00	0.00	1.25	0.50	103.50	5.55	0.80	0.30	02366533
C1125	100.00	1.00	1.25	0.50	104.50	5.55	0.80	0.30	02366541
C1150	100.00	1.00	1.50	0.50	105.00	5.55	0.80	0.30	02366592
C2125	100.00	2.00	1.25	0.50	105.00	5.55	0.80	0.30	02366568
C2150	100.00	2.00	1.50	0.50	106.00	5.55	0.80	0.30	02366606
C3125	100.00	3.00	1.25	0.50	106.50	5.55	0.80	0.30	02366576
C3150	100.00	3.00	1.50	0.50	107.00	5.55	0.80	0.30	02366614
C4125	100.00	4.00	1.25	0.50	107.50	5.55	0.80	0.30	02366584

## MAJOR BENEFITS

- Provides improved anticoagulation properties
- Leads to a higher dose of dialysis
- Demonstrated to reduce heparin need by up to 55%
- Prevents clotting of fibres in the dialyser
- No need to be concerned about hypocalcemia