



Chief Medical Supplies Ltd.

411-19th Street S.E.

Calgary, Alberta T2E 6J7

Telephone: (866) 620-6034

Fax: (888) 620-6074

www.chiefmededical.com

Hemosate Ultra Formulation List

(in mmol/litre after standard 1:45 dilution)

Formulation	Sodium	Potassium	Calcium	Magnesium	Chlorides	Citric Acid	Sodium Acetate	Dextrose	DIN
HS0125	100.0	0.00	1.25	0.5	103.5	0.80	0.30	5.55	02413876
HS1125	100.0	1.00	1.25	0.5	104.5	0.80	0.30	5.55	02413892
HS1150	100.0	1.00	1.50	0.5	105.0	0.80	0.30	5.55	02413930
HS2125	100.0	2.00	1.25	0.5	105.5	0.80	0.30	5.55	02413906
HS2150	100.0	2.00	1.50	0.5	106.0	0.80	0.30	5.55	02413957
HS3125	100.0	3.00	1.25	0.5	106.5	0.80	0.30	5.55	02413914
HS3150	100.0	3.00	1.50	0.5	107.0	0.80	0.30	5.55	02413965
HS4125	100.0	4.00	1.25	0.5	107.5	0.80	0.30	5.55	02413922