



Chief Medical Supplies Ltd.

411-19th Street S.E.

Calgary, Alberta T2E 6J7

Telephone: (866) 620-6034

Fascmile: (888) 620-6074

www.chiefmedical.com

Acid Concentrate Formulation List

(in mmol/litre after standard 1:45 dilution)

Formulation	Sodium	Potassium	Calcium	Magnesium	Chlorides	Acetic Acid	Dextrose	DIN
A1220	100.0	0.00	1.00	0.50	103.0	4.0	8.33	02381257
A1221	100.0	1.00	1.50	0.38	104.8	4.0	5.55	02254344
A1222	100.0	2.00	1.50	0.75	106.5	4.0	5.55	02254336
A1223	100.0	3.00	1.50	0.38	106.8	4.0	5.55	02254328
A1224	100.0	3.00	1.25	0.38	106.3	4.0	5.55	02254301
A1225	100.0	2.00	1.25	0.75	106.0	4.0	5.55	02254298
A1226	100.0	2.00	1.00	0.50	105.0	4.0	5.55	02254255
A1227	100.0	1.00	1.25	0.50	104.5	4.0	5.55	02254247
A1228	100.0	2.00	1.50	0.38	105.8	4.0	5.55	02533774
A1229	100.0	1.00	1.50	0.75	105.5	4.0	5.55	02533901
A1230	104.5	0.00	1.25	0.75	108.5	4.0	8.33	02254085
A1231	104.5	1.50	1.25	0.75	110.0	4.0	5.55	02533855
A1241	102.0	1.00	1.55	0.75	107.6	4.0	5.55	02533782
A1242	102.0	2.00	1.55	0.75	108.6	4.0	5.55	02533847
A1243	102.0	3.00	1.55	0.75	109.6	4.0	5.55	02533758
A1245	102.0	2.00	1.25	0.75	108.0	4.0	5.55	02533812
A1246	102.0	1.00	1.25	0.75	107.0	4.0	5.55	02533820
A1251	104.5	1.50	1.25	0.38	109.3	3.0	5.55	02533790
A1252	104.5	2.00	1.00	0.75	110.0	4.0	5.55	02533839
A1253	104.5	3.00	1.25	0.38	110.8	3.0	11.10	02248823
A1255	104.5	1.50	1.50	0.38	109.8	4.0	5.55	02533766
A1261	100.0	1.00	1.25	0.38	104.3	4.0	5.55	02262479
A1262	100.0	2.00	1.25	0.38	105.3	4.0	5.55	02262487
A1263	100.0	3.00	1.25	0.75	107.0	4.0	5.55	02533871
A1264	100.0	0.00	1.50	0.38	103.8	4.0	11.10	02262509
A1265	100.0	3.00	1.50	0.75	107.5	4.0	5.55	02533863
A1266	100.0	3.00	1.25	0.50	106.5	4.0	8.33	02262525
A1267	100.0	2.00	1.25	0.50	105.5	4.0	5.55	02533898
A1268	100.0	1.00	1.00	0.50	104.0	4.0	8.33	02262541
A1269	100.0	2.00	1.00	0.50	105.0	4.0	8.33	02262568
A1270	100.0	0.00	0.50	0.50	102.0	4.0	5.55	02272512
A1271	110.0	0.00	1.25	0.50	113.5	4.0	5.55	02272520
A1272	100.0	0.00	1.50	0.38	103.8	4.0	5.55	02272539
A1273	100.0	0.00	1.50	0.75	104.5	4.0	5.55	02272547
A1274	100.0	1.00	1.00	0.50	104.0	4.0	5.55	02272555
A1275	100.0	1.50	1.25	0.38	104.8	4.0	5.55	02272563
A1276	100.0	2.00	1.12	0.38	105.0	4.0	5.55	02272571
A1277	100.0	2.00	1.25	0.38	105.3	4.0	5.55	02272598
A1278	110.0	2.00	1.25	0.50	115.5	3.0	5.55	02272601
A1279	100.0	2.00	1.50	0.50	106.0	4.0	5.55	02272628
A1280	100.0	2.50	1.25	0.50	106.0	4.0	5.55	02272636
A1294	100.0	4.00	1.25	0.38	107.3	4.0	5.55	02268728
A1295	100.0	4.00	1.50	0.50	108.0	4.0	5.55	02435888
A1296	100.0	4.00	1.25	0.50	107.5	4.0	5.55	02435918
A1297	100.0	1.00	1.50	0.50	105.0	4.0	5.55	02435969
A1298	100.0	3.00	1.50	0.50	107.0	4.0	5.55	02435950



Chief Medical Supplies Ltd.

411-19th Street S.E.

Calgary, Alberta T2E 6J7

Telephone: (866) 620-6034

Fax: (888) 620-6074

www.chiefmededical.com

Hemosate Ultra Formulation List

(in mmol/litre after standard 1:45 dilution)

Formulation	Sodium	Potassium	Calcium	Magnesium	Chlorides	Citric Acid	Sodium Acetate	Dextrose	DIN
HS0125	100.0	0.00	1.25	0.5	103.5	0.80	0.30	5.55	02413876
HS1125	100.0	1.00	1.25	0.5	104.5	0.80	0.30	5.55	02413892
HS1150	100.0	1.00	1.50	0.5	105.0	0.80	0.30	5.55	02413930
HS2125	100.0	2.00	1.25	0.5	105.5	0.80	0.30	5.55	02413906
HS2150	100.0	2.00	1.50	0.5	106.0	0.80	0.30	5.55	02413957
HS3125	100.0	3.00	1.25	0.5	106.5	0.80	0.30	5.55	02413914
HS3150	100.0	3.00	1.50	0.5	107.0	0.80	0.30	5.55	02413965
HS4125	100.0	4.00	1.25	0.5	107.5	0.80	0.30	5.55	02413922